

Flasks

All laboratories need Flasks of different types. **BOROSIL Scientific** offers a wide range of Flasks for a variety of applications such as boiling, culture, distilling, Erlenmeyer, Filtering, Iodine, and Kjeldahl. Choose one to suit your needs.

BOROSIL Scientific brand Flasks are designed to deliver the best combination of thermal shock resistance, physical strength and resistance to chemical attack. They are made so as to withstand the demands of the modern laboratory. Flasks with interchangeable joints are particularly suitable for all glass apparatus assemblies.

BOROSIL Scientific offers the interchangeable joints most commonly specified for a particular capacity. Tube Adaptors Category Nos. 8800 and 8820 can be used to connect the flask with a particular size joint to any other apparatus having a different size joint.


Made from 3.3 Low Expansion Borosilicate Glass



4060 - FLAT BOTTOM FLASK

Narrow Mouth

- Complies with IS 1381 (Part - I), ISO 1773
- Uniform wall thickness distribution makes these flasks ideal for heating
- All sizes have a durable matt-finished enamelled area for marking



Product Code	Capacity ml	Approx O.D. x Height mm	Approx Neck O.D. mm	Quantity Per Case
4060009	25**	42 x 90	20	50
4060012	50	51 x 90	26	50
4060016	100	64 x 105	26	80
4060018	150**	75 x 132	28	10
4060021	250	85 x 140	34	80
4060022	*300**	88 x 160	24	20
4060024	500	105 x 165	34	40
4060029	1000	131 x 190	42	20
4060030	2000	166 x 230	50	2
4060031	3000	186 x 300	50	4
4060033	5000	224 x 320	58	2
4060034	6000	236 x 315	65	1
4060038	10000	279 x 360	65	1
4060040	20000	350 x 495	78	1


* Suitable for use in Polenske Apparatus

** Not covered in standards

4061 - FLASKS

Flasks, Boiling, Flat Bottom, ISO 24450, Wide Neck

Product Code	Capacity ml	Approx O.D. x Height mm	Approx Neck O.D. mm	Quantity Per Case
4061012D	50	51 x 90	34	50
4061016D	100	64 x 105	34	80
4061021D	250	85 x 140	50	80
4061024D	500	105 x 165	50	40
4061029D	1000	131 x 190	50	20
4061030D	2000	166 x 230	50	2



4100 - FLAT BOTTOM FLASK

Narrow Mouth, Short Neck With I/C Joint

- Uniform wall thickness distribution makes these Flasks ideal for heating
- All sizes have a durable matt-finished enamelled area for marking
- Equipped with interchangeable joint
- The flat base ensures that the Flask can be placed down even without the supporting ring



Product Code	Capacity ml	Interchangeable Joint	Approx Height mm	Quantity Per Case
4100018	150	24 / 29	115	5
4100021	250	24 / 29	130	5
4100024	500	24 / 29	160	5
4100029	1000*	34 / 35	288	2
4100030	2000*	34 / 35	321	2
4100031	3000*	34 / 35	345	2
4100033	5000*	34 / 35	381	2

* These are specific flasks with long neck 34/35 joint used in jewellery manufacturers