



Laboratory Reagent bottles

Precision in every drop

**Made from 3.3 Low Expansion
Borosilicate Glass**



We specialize in providing premium-quality reagent bottles designed to meet the demanding needs of modern laboratories. Our reagent bottle are crafted from high-grade materials such as the borosilicate glass & ensure durability, safety, & long-term chemical stability.


Laboratory Bottle

One essential component found in every laboratory is the **reagent bottle**. From chemicals and solutions to powders and light-sensitive compounds, these sturdy, reliable containers are trusted to safely store a wide range of materials. Though often overlooked, reagent bottles play a **vital role in maintaining accuracy, safety, and efficiency** in daily lab operations. Their robust design, chemical resistance, and versatility make them an **indispensable part of any scientific setup**.

- All bottles from 100 mL onwards feature a standardized screw thread, ensuring the full interchangeability of screw caps & the pouring rings across the range.
- The Bottles, pouring rings and caps are autoclavable and sterilizable.
- For the Amber bottles, ambering is done on the external surface only, hence properties within the bottle remains unchanged.



Industrial applications



1501 & 1519 Reagent Bottle (Clear)


Narrow Mouth With Screw Cap



4796 - 1



1388 (Part 1)



Mechanically strong & Chemical resistant



Ideal for use with light sensitive media and for the long term storage of substances



Scan for Product Details (Clear)



Scan for Product Details (Amber)

Available in the following capacities:

10 ml, 25 ml, 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml, 2000 ml, 3000 ml, 5000 ml, 10000 ml, 20000 ml

1500 & 1509 Reagent Bottle (Clear)

Narrow Mouth with I/C Glass Stopper



Mechanically strong & Chemical resistant



Ideal for use with light sensitive media and for the long term storage of substances



Scan for
Product Details
(Clear)



Scan for
Product Details
(Amber)

Available in the following capacities:
60 ml, 100 ml, 125 ml, 250 ml, 500 ml,
1000 ml, 2000 ml

1506 & 1507 Square Bottle (Clear)

Narrow Mouth with I/C Glass Stopper



Scan for
Product Details
(Clear)



Ideal for Storage & transfer of pharma & biopharma products

Available in the following capacities:
500 ml, 1000 ml, 2000 ml, 5000 ml