

BTI-06 B.O.D. INCUBATOR



- Double-walled construction.
- Outer body made of mild steel powder-coated.
- Inner chamber made of Stainless steel 304.
- **Insulation:** Puff Insulation between two walls.
- **Temp Range:** 5°C to 60°C.deg
- **Temp Accuracy:** +/- 0.5°C.deg
- **Temp controlled by:** PID Micro process controller with dual display to view set & Current temp.
- **Air circulation By:** Air circulation fan motor with a blower to maintain a uniform temperature.
- **Door:** Inner Acrylic door to view sample & outer metal door powder-coated.
- **Gasket:** Silicon Gasket for door.
- **Cooling System:** By a hermetically sealed compressor KIRLOSKER make with CFC-free gas.
- **Wheel:** Castor wheel at the bottom for easy mobility.

Inner working size of the incubator (Capacity)	Liter (Cu.Ft.)	No. of shelves	Internal Dimensions (W x D x H)	Outer Dimensions
3 cu ft	90	2	45 x 45 x 50 cms	24W x 34D x 50H
6 cu ft	170	2	50 x 50 x 70 cms	25W x 38D x 58H
10 cu ft	280	3	50 x 50 x 70 cms	27W x 40D x 66H
12 cu ft	330	3	60 x 60 x 90 cms	29W x 42D x 66H