

Biohazard Bags

- Used for safe handling of infectious and hazardous waste
- Makes waste non-infectious after autoclaving
- Can be autoclaved up to 135°C (275°F) PP only
- For Hazardous Material, autoclave at 121°C (250°F)
- Made of polypropylene, HM.
- Follows Biomedical Waste and Plastic Waste Management Rules
- For high-temperature sterilization, use medical-grade autoclavable bags made from heat-resistant materials like polypropylene

Category: [Lab Utility & Safety](#)



Description

Biohazard Bags

Labsphere presents high-quality **Biohazard Bags** designed for the safe handling, collection, and disposal of infectious and hazardous waste. These bags are manufactured using premium polypropylene material, ensuring durability, reliability, and compliance with biomedical waste management standards.

As a **leading manufacturer of laboratory and cleanroom equipment**, Labsphere provides superior-quality biohazard waste disposal solutions suitable for hospitals, laboratories, research facilities, and pharmaceutical industries. We are a **manufacturer, supplier, dealer, and exporter based in Pune, India**, and we supply across all states in India.

Key Features

- Used for safe handling of infectious and hazardous waste
- Makes waste non-infectious after autoclaving
- Can be autoclaved up to 135°C (275°F) – PP only
- For hazardous material, autoclave at 121°C (250°F)
- Made of high-quality polypropylene (HM)
- Complies with Biomedical Waste and Plastic Waste Management Rules
- Suitable for high-temperature sterilization
- Strong, leak-proof, and durable construction

Applications

- Hospitals and healthcare facilities
- Microbiology and pathology laboratories
- Pharmaceutical industries
- Cleanroom environments
- Research and diagnostic centers
- Biomedical waste disposal units

Technical Details

- Material: Polypropylene (HM)
- Temperature Resistance: Up to 135°C (PP)
- Autoclaving Standard: 121°C for hazardous material
- Compliance: Biomedical Waste Management Rules
- Available Colors: Red, Yellow, Black

Code	Sizes (WxL inch)	Color	Pcs per/pack	HSN No.
LS25-5001	8x12	Red	200	39.23.2990
LS25-5002	12x24	Red	200	39.23.2990
LS25-5003	12x18	Red	200	39.23.2990
LS25-5004	18x24	Red	200	39.23.2990
LS25-5005	25x35	Red	100	39.23.2990
LS25-5006	42x42 (HM)	Red	50	39.23.2990
LS25-5007	37x48 (HM)	Red	50	39.23.2990
LS25-5152	28x43 (HM)	Red	50	39.23.2990
LS25-5008	8x12	Yellow	200	39.23.2990
LS25-5009	12x24	Yellow	200	39.23.2990
LS25-5010	12x18	Yellow	200	39.23.2990
LS25-5011	18x24	Yellow	200	39.23.2990
LS25-5012	25x35	Yellow	100	39.23.2990
LS25-5013	42x42 (HM)	Yellow	50	39.23.2990
LS25-5014	42x42 (HM)	Black	50	39.23.2990
LS25-5015	37x48 (HM)	Yellow	50	39.23.2990
LS25-5151	28x43 (HM)	Yellow	50	39.23.2990