

NEW**5**
Year
Warranty

Digital Vortexer

Controlled, high-speed touch mixing

The Digital Vortexer combines digital precision, ergonomic comfort, and durable construction, featuring variable speed control, low activation force, high-speed mixing, stability, chemical resistance, compliance certifications and a 5-year warranty

Digital Control: produces more accurate results

- Digital controls enhance user control for precise mixing
- Illuminated digital display ensures accuracy, reproducibility, and consistent results
- Variable speed control offers flexibility and precise control of mixing intensity
- Ergonomic: prevents user fatigue
- Minimal activation force, supporting both touch start and continuous operation for convenience
- High-speed vortex capability – up to 3,600 RPM
- Multiple-position head allows vortexing at angles beyond 90 degrees for enhanced versatility
- Durable Design: guarantees quality that saves money
- Constructed from ABS and silicone, unit resists common chemicals and facilitates easy cleaning
- Stable platform and elastomer feet ensure unit stays stable during operation
- Meets Safety and Compliance Certifications while offering lower energy consumption
- 5-Year warranty



Item No.	Color	L x W x H in	L x W x H cm
120805	Grey/Purple	7.9 x 6.1 x 5.6	20.1 x 15.5 x 14.2
120806	Grey/Blue		

Specifications

Speed Range	500 – 6,000 rpm
Operating Modes	Intermittent (Touch) and Continuous
Intermittent Speed Range	600 to 3,600 RPM
Continuous Speed Range	600 to 2,600 RPM
Run Time Range	0 to 50 hrs
Timer Range:	Down-counting, in minutes and seconds, (Continuous Mode run time is up to 50 hours, Intermittent Mode run time is up to 1 minute), Upward-counting run time is up to 50 hrs
Weight	5.4 lbs (2.4 kg)
Electrical Input	100-240 VAC, 2.5A, 50/60 Hz Output: 12 VDC, 3.0A, 36.0W
Maximum Load (Continuous)	250 grams (0.55 lbs)