



e count

Rapid, accurate colony counter

Combines a multi-function electronic counter with a Sharpie® Pen for marking to prevent missing or double-counting colonies

- Light, ergonomic body is balanced for easy manipulation
- Select between an audible beep or an LED flash to acknowledge count
- Store up to 32 separate counts and display the total of all counts
- Count up/down verification
- Quick and easy to change between different colored markers to check different colony counts in the same view

Package Includes:

- eCount Colony Counter
- Black Sharpie® Pen (Page 35)
- Long-life button battery included
- Hand tally Adapter
- Weighted Stand



US Patent No. D564,385

EP Registered Design No. 705751-0001

Item No.	L x W in	L x W cm
120000	7.4 x 1.2*	17 x 3*
120228	CR1225 Battery	

*Measurements shown are with Sharpie® Pen inserted



eCount with Mini-Light Box
with Mini-Magnifier, page 31



Weighted Stand
for bench top use



Battery and Pen not
included in warranty



Lock Counter Off

to allow use of the marker without
registering more counts

Hand Tally Adapter
allows you to use the
eCount as a digital tally
counter

