

BTI-25 INCUBATORS BACTERIOLOGICAL LAB TYPE



- Construction: Double-walled.
- Outer body made of mild steel powder-coated.
- Inner chamber made of stainless steel 304 grade.
- **Insulation:** Glass wool insulation between two walls.
- **Temp range:** 37°C to 60 °C (To maintain 37°C, the unit should be operated in an environment with a room Temperature of 30°C or below)
- **Accuracy:** $\pm 0.5^{\circ}\text{C}$.
- Temp Controlled by: Microprocessor PID controller with dual display & auto tune facility.
- With G.I. wire Mesh shelves
- **Door:** 2-door inner acrylic door to view sample & outer metal door.
- **Heater:** special Nichrome heater designed for long life & saving of power consumption.

Inner working size in Inches	Outer Dimension in inches	Rating Watts	Capacity	No of shelves
12"x12"x12"	18W x 19D x 28H	250	27 LTR	1
14"x14"x14"	21W x 21D x 28H	300	44 LTR	2
18"x18"x18"	25W x 25D x 33H	400	90 LTR	2