

DURAN® Original GL 25 Laboratory Bottle, clear, with screw cap, 10 mL

- Complies with the standards in ISO 4796-1:2016 Laboratory glassware -- Bottles -- Part 1: Screw-neck bottles (does not apply to 10 ml bottle)
- Manufactured from DURAN® borosilicate glass 3.3. DURAN® is a neutral glass, classified as a Glass Type I according to USP <660>, EP (3.2.1) and JP <7.01>, and meets ASTM E438, Type I, Class A requirements
- Standardised GL-thread and matching cap system for particularly tight sealing and simple, clean pouring out
- Retrace Code, and certificate for manufacturing lot traceability download at <https://cert.dwk.com/>
- Fired-on, highly durable white print
- Suitable for long-term storage
- Easy-to-read scale with graduation marks to estimate the content

Catalog Number: **218010851**

Product Overview

DURAN® Original GL laboratory bottles have become the gold standard of multifunctional laboratory bottles since their introduction in 1972. The bottles are moulded from high-purity, clear, USP/EP Type I borosilicate glass 3.3, renowned for its excellent thermal performance, and high resistance to chemical attack. A combination of high-performance materials and robust ISO 4796-1 design provides a long working life. The bottles feature an easy-to-read graduated volume scale, and marking field in permanent white enamel. Permanent Retrace Code on each bottle for manufacturing lot traceability. DURAN® Original GL bottles are supplied complete with iconic blue polypropylene cap and drip-free pouring ring for use up to 140 °C. All components are fully autoclavable at 121 or 134 °C.

Specifications

Boro 3.3

Boro 3.3 [See Details](#)

Remarks	Does not conform with ISO 4796-1 standard. With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).
---------	--

Material	Boro 3.3, PP
----------	--------------

Graduated	Yes
-----------	-----

Color	Clear
-------	-------

Colour	Clear
--------	-------

Diameter d (mm)	33
-----------------	----

Height h (mm)	59
---------------	----

Capacity (mL)	10
---------------	----

Brim Capacity (mL)	20
--------------------	----

Diameter (mm)	33
---------------	----

Thread	GL 25
--------	-------

Glass Color	Clear
-------------	-------

Marking Area	Marking Area
--------------	--------------

Height (mm)	59
-------------	----

Autoclavable	yes
--------------	-----

Retrace Code	Retrace Code
--------------	--------------

Quantity	10
----------	----

Weight (kg)	0.037
-------------	-------