

BTI-31 LAMINAR AIR FLOW CABINETS



- **Construction:** Outer body made of M.S powder coated & Inner working area made of S.S.
- **HEPA filter:** Efficiency of HEPA filter of 99.97 down on 0.3 microns & pre filter down to 5 microns.
- A suitable blower fitted at the back of the uniform flow.
- Side Acrylic panel.
- Door fitted with an acrylic door.
- Supplied with 1 U.V light & 1 day light.
- Supplied with 1 analogue manometer to measure pressure flow.
- Supplied with a gas cork.

HORIZONTAL TYPE

Model type	Inner working Size in feet	Outer dimension in inches
HM-1	2" x 2 x 2"	28"(W) x 37"(D) x 57"(H)
HM-3	3" x 2" x 2"	39"(W) x 39"(D) x 57"(H)
HM-4	4" x 2" x 2"	52"(W) x 42"(D) x 57"(H)

VERTICAL TYPE

Model type	Inner working Size in feet	Outer dimension in inches
VM-2	2" x 2 x 2"	28"(W) X 28"(D) x 80"(H)
VM-3	3" x 2" x 2"	39"(W) x 28"(D) x 80"(H)
VM-4	4" x 2" x 2"	52"(W) x 28"(D) x 80"(H)