

BTI-36 MUFFLE FURNACE



- Double-walled construction.
- Outer body made of Mild steel.
- Inner chamber made of ceramic Muffle KUMAR Make.
- **Insulation:** A special type of ceramic Insulation between the body.
- **Coil:** Kanthal A1 Swedish-made coil for long-lasting.
- **Temp Range:** 400 C to max 1150 °C.
- **Accuracy:** +/- 3 deg.
- Temp Controlled by: PID Dual display to view set & current temp.
- **Safety:** A special type of thermocouple fuses to protect against overheating.
- **Sensor:** 'K' type calibrated sensor.

Inner working Size in inches	Outer dimension in inches	Rating watts
9"x4"x4"	22"D x 15"Wx 23"H	2000
10"x5"x5"	24"D X 16" WX 24"H	2500
12"x6"x6"	25D"X 17"WX 26"H	3500

SPARE ACCESSORIES

Particulars
Thermo-Couple
Electronic Digital Timer

SPARE MUFFLE WITH COIL

Size in inches	Rating
9" x 4" x 4"	2.0 KW
10" x 5" x 5"	2.5 KW
12" x 6" x 6"	3.5 KW