

Nylon Membrane Filter

Nylon membrane filters for HPLC and lab filtration, suitable for aqueous and organic solutions. Available in multiple sizes and pore options. Manufactured by Labsphere, Pune, India.

Category: [Filtration](#)



Description

Nylon Membrane Filter

The **Nylon Membrane Filter** is a high-performance filtration solution designed for effective particle removal in laboratory, cleanroom, and pharmaceutical applications. These filters are ideal for aqueous solutions and organic solvents, making them suitable for analytical processes such as HPLC and routine laboratory filtration.

Manufactured by **Labsphere**, we are a leading manufacturer of laboratory and cleanroom equipment, and a trusted manufacturer, supplier, dealer, and exporter based in Pune, India. We supply across all states in India, delivering premium-quality filtration products and microbiology lab supplies.

Key Features

- Recommended for particle removal in water, aqueous solutions, and organic solvents
- Ideal for HPLC and analytical filtration applications
- High chemical compatibility and mechanical strength
- Available in multiple sizes for versatile use
- Ensures consistent and reliable filtration performance

Applications

- HPLC sample preparation
- Laboratory filtration processes
- Pharmaceutical and research applications
- Microbiology and analytical laboratories
- Chemical and environmental testing

Technical Details

- Membrane material: Nylon
- Available sizes: 13mm, 25mm, 47mm (additional sizes 90mm, 142mm, 293mm available)
- Pore sizes: 0.22µm, 0.45µm (also available in 0.8µm, 1.2µm, and 5.0µm)
- High flow rate and efficient filtration

Code	Size	Pcs per/pack	HSN No.
LS23-4079	13mm, 0.22µm	100	84.21.0000
LS23-4080	25mm, 0.22µm	100	84.21.0000
LS23-4081	47mm, 0.22µm	100	84.21.0000
LS23-4082	13mm, 0.45µm	100	84.21.0000
LS23-4083	25mm, 0.45µm	100	84.21.0000
LS23-4084	47mm, 0.45µm	100	84.21.0000