

1500 - REAGENT BOTTLES


Plain Narrow Mouth, Graduated with Inter-changeable Flat Head Stopper

- Complies with IS 1388 (Part II), EN ISO 4796 - 2
- These narrow mouth reagent storage bottles have a sturdy ground neck
- All sizes have a durable matt-finished enamelled area for marking
- These bottles are mechanically strong and chemical resistant
- Supplied with solid glass stopper

Product Code	Capacity ml	Approx O.D. x Height mm	Neck Stopper Size	Quantity Per Case
1500013	60**	45 x 100	14 / 23	40
1500016	100	56 x 105	14 / 23	10
1500017	125*	56 x 145	19 / 26	50
1500021	250	70 x 145	19 / 26	10
1500024	500	86 x 180	24 / 29	10
1500029	1000	100 x 230	29 / 32	10
1500030	2000	136 x 255	34 / 35	10

- *Height indicated is of Bottle only

- ** 60 & 125 ml is not covered in standards



1501 - REAGENT BOTTLES

Graduated With Screw Cap and Pouring Ring

- Complies with IS 1388 (Part I), EN ISO 4796 - 1
- These bottles are mechanically strong and chemical resistant
- Designed with plastic pouring ring for drip-free operation
- The screw caps and pouring ring are made of PP
- The bottles, screw caps and pouring rings can be sterilized and autoclavable
- Printed with a Retrace Code
- 100 ml to 20000 ml Bottles have uniform GL 45 thread

Product Code	Capacity ml	Approx O.D. x Height* mm	Thread Specification	Quantity Per Case
1501006	10**#	36 X 50	GL 25	10
1501009	25#	36 x 70	GL 25	20
1501012	50	46 x 87	GL 32	20
1501016	100	56 x 100	GL 45	10
1501018	150**	60 x 120	GL 45	10
1501021	250	70 x 140	GL 45	10
1501024	500	86 x 176	GL 45	10
1501029	1000	101 x 225	GL 45	10
1501030	2000	136 x 260	GL 45	10
1501031	3000**	158 x 295	GL 45	2
1501033	5000***	183 x 323	GL 45	1
1501038	10000**	235 x 405	GL 45	1
1501040	20000***	296 x 500	GL 45	1

- *Height indicated is of Bottle only

- # Without Pouring Ring

- ** Not covered in any standard

- *** As per Borosil Dimension