

Reservoirs

Buy high-quality **167 – Laboratory Reservoirs** available in 10 ml, 25 ml, and 50 ml capacities. Designed with graduated walls, angled bottom, and spill-proof spout. Available in bulk, sterile, and individually wrapped options. Trusted manufacturer, exporter, dealer, and supplier in Pune, India with all India delivery.

Category: [Lab Equipments](#)



Description

167 – Reservoirs – Laboratory Liquid Handling Reservoirs

The **167 – Reservoirs** are high-quality laboratory consumables designed for efficient liquid handling and reagent storage. These reservoirs are manufactured from durable **polystyrene**, ensuring stability, clarity, and reliable performance in laboratory environments.

We are a leading **manufacturer, exporter, dealer, and supplier of Laboratory Reservoirs in Pune, India**, supplying premium lab consumables across major cities including Mumbai, Delhi, Bangalore, Hyderabad, Chennai, Kolkata, Ahmedabad, Jaipur, Surat, Lucknow, Patna, Bhopal, Indore, Nagpur, and all states across India.

Key Features of Reservoirs

- **Graduated walls** for quick and accurate liquid measurement
- **Compact polystyrene body** with broad edges for stability
- **Angled bottom design** ensures maximum reagent recovery
- **Pour-off spouts** minimize spillage during use
- **Extra-large stable base** prevents tipping
- Available in **10 ml, 25 ml, and 50 ml capacities**
- Available in **bulk, sterile, and individually wrapped sterile options**

Applications

These reservoirs are widely used in laboratories, research centers, biotechnology labs, pharmaceutical industries, and diagnostic labs for liquid dispensing, reagent handling, and pipetting applications.

Code	Product	Pack Size	HSN No.
LS21-2964	10ml Bulk	100x5	90.31.9000
LS21-2965	10ml Sterile	5x40	90.31.9000
LS21-2966	10ml individually Wrap Sterile	1x50	90.31.9000
LS21-2967	25ml Bulk	100x5	90.31.9000
LS21-2968	25ml Sterile	40x5	90.31.9000
LS21-2969	25ml individually Wrap Sterile	1x50	90.31.9000
LS21-2970	50ml Bulk	500	90.31.9000
LS21-2971	50ml Sterile	40x5	90.31.9000
LS21-2972	50ml individually Wrap Sterile	1x 50	90.31.9000