

6402-FUNNELS

Separating, Pear Shape, Fitted with Boroflo Stopcock with PTFE Key and Interchangeable Stopper

- Complies with IS 1575, DIN / ISO 4800
- Made 3.3 borosilicate glass, USP Type I
- Equipped with a Boroflo stopcock with PTFE key
- Pear shape, with a stem
- Generally used in the separation phase



Product Code	Capacity ml	Stopper Size	Quantity Per Case
6402013	60*	14 / 23	10
6402017	125*	19 / 26	10
6402021	250	24 / 29	10
6402024	500	24 / 29	10
6402029	1000	29 / 32	10
6402030	2000	29 / 32	10


* Not covered in standards.



6403-FUNNELS

Separating, Pear Shape, with PTFE Key Stopcock, Interchangeable Stopper

- Complies with IS 1575, DIN / ISO 4800
- Made of 3.3 borosilicate glass
- USP Type I glass
- Equipped with a PTFE key stopcock
- Pear shape, with a stem
- Generally used in the separation phase



Product Code	Capacity ml	Stopper Size	Quantity Per Case
6403017	125*	19 / 26	10
6403021	250	24 / 29	10
6403024	500	24 / 29	10
6403029	1000	29 / 32	10
6403030	2000	29 / 32	10


* Not covered in standards.



6405 - FUNNEL

Pressure Equilising Cylindrical, with PTFE Stopcock

- Made 3.3 borosilicate glass, USP Type I
- Equipped with totally grease free and self lubricating PTFE stopcock
- Supplied with polypropylene stopper



Product Code	Capacity ml	Socket	Cone	Quantity Per Case
6405012	50	19 / 26	19 / 26	5
6405016	100	19 / 26	19 / 26	5
6405021	250	24 / 29	24 / 29	5
6405024	500	24 / 29	24 / 29	5

