

Safety Bottle Carriers

Secure and convenient bottle transport

- Ergonomic handle prevents hand fatigue and related accidents when carrying heavier loads
- Handle swings down to either side for easier loading and unloading
- Centering ring stabilizes bottle in the center of carrier for safer, steadier transport
- Carriers come in four colors with designated labeling areas for easy identification
- Made of chemical resistant thermoplastic elastomer - PP blend



Safety Bottle Carrier, 4.5 L

Item No.	Color	Dimensions (Dia. x H) in	Dimensions (Dia. x H) cm
120911	Red	7.2 x 10.5 (16 in height, handle extended)	18.3 x 26.7 (40.6 cm height, handle extended)
120912	Black		
120913	Blue		
120914	Green		

Available in 4 Colors
for easy identification and organization

Safety Bottle Carrier, 1.5 L

Item No.	Color	Dimensions (Dia. x H) in	Dimensions (Dia. x H) cm
120915	Red	4.5 x 6.5 (10 in height, handle extended)	11.4 x 16.5 (25.4 cm height, handle extended)
120916	Black		
120917	Blue		
120918	Green		



Movable Handle
swings down to allow
clear access to vessel

