

Syringe Filter-Non sterile Nylon/PVDF/ PTFE/ PES/ CA

Sterile syringe filters available in PVDF, PTFE, and Nylon with 0.22µm and 0.45µm pore sizes. Ideal for lab and pharma filtration. Manufactured by Labsphere, Pune, India.

Category: [Filtration](#)



Description

Syringe Filter – Sterile PVDF / PTFE / Nylon

The **Syringe Filter – Sterile PVDF / PTFE / Nylon** is designed for efficient and reliable filtration in laboratory, cleanroom, and pharmaceutical applications. These high-quality sterile syringe filters ensure precise filtration performance, excellent chemical compatibility, and consistent results for both aqueous and organic solutions.

Manufactured by **Labsphere**, we are a leading manufacturer of laboratory and cleanroom equipment, and a trusted manufacturer, supplier, dealer, and exporter based in Pune, India. We supply across all states in India, delivering premium filtration products and microbiology lab supplies.

Key Features

- Available in PVDF, PTFE, and Nylon membrane types
- High filtration efficiency with 0.22µm and 0.45µm pore sizes
- Compatible with aqueous and organic solutions
- Ensures reliable and consistent filtration performance
- Suitable for laboratory, pharmaceutical, and research applications
- Available in multiple diameters: 13mm, 25mm, and 33mm

Applications

- Sample preparation and clarification
- Laboratory filtration processes
- Pharmaceutical and biotech applications
- Microbiology and analytical labs
- Chemical and research laboratories

Technical Details

- Sterile syringe filter units
- Membrane options: PVDF, PTFE, Nylon
- Pore sizes: 0.22µm and 0.45µm
- Multiple size configurations for versatile use
- Also available in PES / CA variants

Syringe Filter - Sterile PVDF/PTFE/Nylon			
Code	Product	Pcs per/pack	HSN No.
LS23-4031	PVDF 0.22µm, 13mm	50	84.21.0000
LS23-4032	PVDF 0.22µm, 25mm	50	84.21.0000
LS23-4033	PVDF 0.22µm, 33mm	50	84.21.0000
LS23-4034	PVDF 0.45µm, 13mm	50	84.21.0000
LS23-4035	PVDF 0.45µm, 25mm	50	84.21.0000
LS23-4036	PVDF 0.45µm, 33mm	50	84.21.0000
LS23-4037	PTFE 0.22µm, 13mm	50	84.21.0000
LS23-4038	PTFE 0.22µm, 25mm	50	84.21.0000
LS23-4039	PTFE 0.22µm, 33mm	50	84.21.0000
LS23-4040	PTFE 0.45µm, 13mm	50	84.21.0000
LS23-4041	PTFE 0.45µm, 25mm	50	84.21.0000
LS23-4042	PTFE 0.45µm, 33mm	50	84.21.0000
LS23-4043	Nylon 0.22µm, 13mm	50	84.21.0000
LS23-4044	Nylon 0.22µm, 25mm	50	84.21.0000
LS23-4045	Nylon 0.22µm, 33mm	50	84.21.0000
LS23-4046	Nylon 0.45µm, 13mm	50	84.21.0000
LS23-4047	Nylon 0.45µm,25mm	50	84.21.0000
LS23-4048	Nylon 0.45µm, 33mm	50	84.21.0000

Available in: PES/CA

Code	Product	Pcs per/pack	HSN No.
LS23-4001	PVDF 0.22µm, 13mm	100	84.21.0000
LS23-4002	PVDF 0.22µm, 25mm	100	84.21.0000
LS23-4003	PVDF 0.22µm, 33mm	100	84.21.0000
LS23-4004	PVDF 0.45µm, 13mm	100	84.21.0000
LS23-4005	PVDF 0.45µm, 25mm	100	84.21.0000
LS23-4006	PVDF 0.45µm, 33mm	100	84.21.0000
LS23-4007	PTFE 0.22µm, 13mm	100	84.21.0000
LS23-4008	PTFE 0.22µm, 25mm	100	84.21.0000
LS23-4009	PTFE 0.22µm, 33mm	100	84.21.0000
LS23-4010	PTFE 0.45µm, 13mm	100	84.21.0000
LS23-4011	PTFE 0.45µm, 25mm	100	84.21.0000
LS23-4012	PTFE 0.45µm, 33mm	100	84.21.0000
LS23-4013	Nylon 0.22µm, 13mm	100	84.21.0000
LS23-4014	Nylon 0.22µm, 25mm	100	84.21.0000
LS23-4015	Nylon 0.22µm, 33mm	100	84.21.0000
LS23-4016	Nylon 0.45µm, 13mm	100	84.21.0000
LS23-4017	Nylon 0.45µm,25mm	100	84.21.0000
LS23-4018	Nylon 0.45µm, 33mm	100	84.21.0000
LS23-4019	PES 0.22µm, 13mm	100	84.21.0000
LS23-4020	PES 0.22µm, 25mm	100	84.21.0000
LS23-4021	PES 0.22µm, 33mm	100	84.21.0000
LS23-4022	PES 0.45µm, 13mm	100	84.21.0000
LS23-4023	PES 0.45µm, 25mm	100	84.21.0000
LS23-4024	PES 0.45µm, 33mm	100	84.21.0000
LS23-4025	CA 0.22µm, 13mm	100	84.21.0000
LS23-4026	CA 0.22µm, 25mm	100	84.21.0000
LS23-4027	CA 0.22µm, 33mm	100	84.21.0000
LS23-4028	CA 0.45µm, 13mm	100	84.21.0000
LS23-4029	CA 0.45µm, 25mm	100	84.21.0000
LS23-4030	CA 0.45µm, 33mm	100	84.21.0000