

BTI – 53 (A) VACCUM PUMP FOR VACCUM OVEN



- Model No: 550-Z
- Max. Flow (Ltrs/Min): 45
- Max. Vacuum (inch): 29"
- Max. Press (PSI_g): 55
- Approx. Weight (kg): 9.0
- Motor HP: 3/4
- Approx. Dimension (mm): 290 x 120 x 180
- Note: The warranty of 1 year on the vacuum pumps/compressor is only against any manufacturing defects. No warranty if foreign material enters the pump through the suction port and damages the pump.

BTI-53 (B) VACCUM PUMP (OIL FREE DIAPHRAM TYPE)



TID 15



TID 25



TID 45



TID 75 S

MODEL	Max Flow (Ltrs/Min)	Max Vacuum (inches Hg)	Max Press (PSI _g)	Motor HP
TID-15	15	22"	25	1/20
TID-25-S	15	27"	35	1/16
TID-45	45	22"	35	1/8
TID-75-S	45	27"	55	1/4

- NOTE: THE WARRANTY OF 1 YEAR ON THE VACUUM PUMPS/COMPRESSOR IS ONLY AGAINST ANY MANUFACTURING DEFECTS. NO WARRANTY IF FOREIGN MATERIAL ENTERS THE PUMP THROUGH THE SUCTION PORT AND DAMAGES THE PUMP.