



## 7102 - PIPETTES

### Transfer, Volumetric, Class B

- Complies with EN ISO 648
- Calibrated for delivery
- Generally used for accurate measurement and decanting of liquid

| Product Code | Capacity ml | Tolerance $\pm$ ml | Quantity per Case |
|--------------|-------------|--------------------|-------------------|
| 7102001D     | 1           | 0.015              | 10                |
| 7102002D     | 2           | 0.02               | 10                |
| 7102005D     | 5           | 0.03               | 10                |
| 7102006D     | 10          | 0.04               | 10                |
| 7102008D     | 20          | 0.06               | 10                |
| 7102009D     | 25          | 0.06               | 10                |
| 7102012D     | 50          | 0.10               | 10                |
| 7102016D     | 100         | 0.15               | 10                |



## 7103 - PIPETTES

### Transfer, Volumetric, Class B

- Complies with ASTM E-969-02
- Calibrated for delivery
- Generally used for accurate measurement and decanting of liquid

| Product Code | Capacity ml | Tolerance $\pm$ ml | Quantity per Case |
|--------------|-------------|--------------------|-------------------|
| 7103001      | 1           | 0.012              | 10                |
| 7103002      | 2           | 0.012              | 10                |
| 7103003      | 3           | 0.02               | 10                |
| 7103004      | 4           | 0.02               | 10                |
| 7103005      | 5           | 0.02               | 10                |
| 7103706      | 6           | 0.03               | 10                |
| 7103707      | 7           | 0.03               | 10                |
| 7103708      | 8           | 0.04               | 10                |
| 7103709      | 9           | 0.04               | 10                |
| 7103006      | 10          | 0.04               | 10                |
| 7103007      | 15          | 0.06               | 10                |
| 7103008      | 20          | 0.06               | 10                |
| 7103009      | 25          | 0.06               | 10                |
| 7103010      | 30*         | 0.06               | 10                |
| 7103011      | 40*         | 0.10               | 10                |
| 7103012      | 50          | 0.10               | 10                |
| 7103016      | 100         | 0.16               | 10                |

\* Not covered in ASTM standards