

## BTI-57 WATER BATH THERMOSTATIC - RECTANGULAR



- Double-walled construction.
- Outer body made of stainless steel grade 304.
- Inner body made of stainless steel 304 grade.
- **Temperature range:** 40°C to 99°C.
- Temperature controlled by a hydraulic thermostat.
- **Temperature accuracy:**  $\pm 3^\circ\text{C}$  to  $4^\circ\text{C}$ .
- Supplied with a 7.5 cm concentric ring.
- Fitted with a concentric ring cover on top.
- Fitted with a concentric ring cover on top of the lid to keep a 250 or 500 ml flask on the lid.
- **Note:** The heating element is not covered under warranty against failure caused by dry run operation or low water level.

| Inner working size<br>in CM | Rating watts | No. of Holes | Capacity in LTR |
|-----------------------------|--------------|--------------|-----------------|
| 35 x 25 x 10 cm             | 1000         | 6 holes      | 8 LTR           |
| 40 x 30 x 10 cm             | 1500         | 12 holes     | 12 LTR          |